

# Medical Gauze Dressing

## 一、Product name, model and specification

Product name: Medical gauze dressing

Model: type A, gauze block; Type B, abdominal towel (excluding abdominal towel band); Type C, gauze cotton pad; Type D, gauze wrap. A, B, C, D also have conventional and x-ray detectable type, sterilization and non-sterilization supply.

Specification: See the product technical requirements for details

## 二、Name, address, production address, contact information and after-sales service unit of the registered person/manufacturer

- 1、Name: Hunan Zhenheyikang Medical Instrument Co.,LTD
- 2、Address: No.6 Building Jingxiang Energy, No.55 Xianguang east road, Gaoxin district Xiangtan, Hunan
- 3、Production address: NO.6 Building, NO.8 building, 3<sup>rd</sup> and 4<sup>th</sup> floor Jingxiang Energy, No.55 Xianguang east road, Gaoxin district Xiangtan, Hunan, NO 19 Chutian Road, Yuetang District, Xiangtan City, Hunan
- 4、Contact number: +86 17711627084
- 5、After-sales service unit: Hunan Zhenheyikang Medical Instrument Co.,LTD

三、Manufacturing License NO: Hunan Medical Products Administration Production Licence No. 20150013

四、Product Registration Certificate No, Product Technical Requirements No: HNFDA Certified 20182140093

## 五、Structure and composition, scope of application

The product adopts skim cotton gauze, skim cotton adhesive mix spinning gauze cloth, x-ray detectable gauze as raw materials, can be dyed, the product is made of direct gauze material or cotton by folding or sewing.

It is used to clean the skin, mucous membrane or wound surface, and to protect the wound surface together with drugs commonly used in wound care. It can also be used to absorb internal exudate, compress hemostasis and support, and protect organs and tissues during surgery.

六、Use method: open the bag can be used

## 七、Contraindications, precautions, warnings and reminders

Product contraindications: None

Damaged packaging prohibited use, expired products prohibited use;

Disposable use, please destroy after use

Avoid high temperature, humidity, direct sunlight

## 八、Interpretation of graphics, symbols, abbreviations and other contents used in medical device labels




一次性使用



参考使用说明书



包装破损切勿使用



环氧乙烷灭菌



无菌

## 九、Installation and use instructions or illustrations

This product does not need to be installed, open the bag can be used.

## 十、Product maintenance and maintenance methods, special storage, transportation conditions, methods

Products should be stored in relative humidity  $\leq 80\%$ , non-corrosive gas and well-ventilated indoor, ordinary transport, in the transportation, storage process should be kept clean and pollution-free, strictly prohibited rain and snow immersion.

## 十一、Date of production, shelf life or expiration date

Valid for three years, production batch number, production date, expiry date please refer to the small package or qualification label.

**十二、 Accessories list, including accessories, appendix, lossitage replacement cycle and replacement method description**

This product has n o accessories

**十三、 Date and version of preparation or revision of the specification**

This specification was revised in February 2023, version 2023.